





Beginner's Manual of Experiments at SACLA


September 10, 2019

Written by Y. Kubota (kubota@spring8.or.jp)

Summary of data handling at SACLA

Run#	Tag#	Simple data <ul style="list-style-type: none"> pulse motor (PM) position digitizer readout (PD, BM, etc) preamp gain shutter status etc 	Image data <ul style="list-style-type: none"> MPCCD IMPERX OPAL
800051	1771706001	PM PD shutter . . .	
	1771706002	PM PD shutter . . .	
	⋮	⋮	⋮
(*)	1771706011	PM PD shutter . . .	not available
	1771706012	PM PD shutter . . .	
	1771706013	PM PD shutter . . .	
800052	1771706014	PM PD shutter . . .	
	1771706015	PM PD shutter . . .	
	⋮	⋮	⋮

(*) changing conditions such as motor position

automatic

DB

on-demand

Storage

◆ What is “Tag” and “Run”?

- **Tag** . . . Most data (**except for image data**) are recorded **automatically** on Database (DB) shot-by-shot. Each of the shots is identified by a unique Tag number.

- **Run** · · · An index to handle a lump of shots. Each run is identified by a unique Run number.

- ◆ Use “RunControl”, when you want
 - To record image.
 - To make Runs.
 - To scan with more than 1 axis or to execute complex scans.

- ◆ To extract data from DB and Storage, refer to the Run numbers to obtain a list of Tag numbers.

Summary of how to make a scan

- ① **Configure scan parameters**
 - Use “ScanConfigEditor” GUI.
 - Write the scan conditions into a csv file.
- ② **Set up RunControl**
 - Load the csv file into RunControl
 - Select the image detectors.
 - Input the number of XFEL shots.
- ③ **Start scan in RunControl**
- ④ **Extract data from the DB and the Storage**
 - Use “DataConverterGUI”.
 - Extract data by referring the Run numbers allocated in ③.
- ⑤ **Analyze the data**

Good luck with analysis.